





Standard Operating Procedure ANTHRAX



PARTNER Christian Medical College (CMC), Vellore and National Center for Disease Control (NCDC) Ministry of Health and Family Welfare

Acknowledgements:

Dr BalajiVeeraraghavan, Professor, Department of Microbiology, CMC, Vellore for his support in preparing this document.

The Technical Advisory Committee for Anthrax initiative for their guidance.

National Referral Laboratories for Diagnosis of Human Anthrax in India

- 1. The Department of Microbiology, Christian Medical College (CMC), Vellore
- 2. Centre for Arboviral and Zoonotic diseases, The National Center for Disease Control (NCDC), Ministry of Health and Family Welfare



© Indian Council of Medical Research 2020. All rights reserved Publications of the ICMR can be obtained from ICMR-HQ, Indian Council of Medical Research Division of Epidemiology and Communicable Diseases (ECD), Ansari Nagar, New Delhi, India

Request for permission to reproduce or translate ICMR publications whether for sale ornoncommercial distribution should be addressed to ICMR-Headquarters at the above address. The mention of specific companies or certain manufacturers' products does not imply that theyare endorsed or recommended by the ICMR in preference to others of a similar nature that arenot mentioned. The responsibility for the interpretation and use of the material lies with the reader. In no eventshall ICMR be liable for damages arising from its use. For any suggestions, please write to us at ICMR.





प्रोफेसर (डा.) खलराम भागेव, एटव श्र एवएस. क्रांस एक्सास्वेव (से) एक्सासोव से) एवएसते एवएस. एवएस्स एक्सास्य अनुसंधन विभाग साचिव, भारत सरकार स्वाप्य व्य परियाद करवाग मकारण एवं महानिदेशक, आई सी एम आर Prof. (Dr.) Balram Bhargava, Padma Shri M.O. DM. FRCP (Glass.). FRPC (Edn.). FACC, FAMA, FMAS, FMAS, FASc

Secretary to the Government of India Department of Health Research Ministry of Health & Family Welfare & Director-General, ICMR





भारतीय आयुर्विज्ञान अनुसंधान परिषद स्वाख्य अनुसंधान विमाग स्वाख्य एवं पश्चिर कल्याण मंत्रात्मय वै, वमविंगस्वामी भवन, अंसरी नगर नई दिल्ली - 110 029 (मारव) Indian Council of Medical Research Department of Health Research Ministry of Health & Family Welfare V. Ramalingaswami Bhawan, Ansari Nagar New Delni - 110 029 (INDIA)

FOREWORD

Anthrax is a neglected disease with a potential for natural or manmade outbreaks. The extraordinary stability of B.anthracis spores that survive in soil for around 200 years and their resistance to many of the disinfection procedures demands extremely careful handling of samples or bacterial cultures and following stringent decontamination procedures.

Anthrax remains a remarkably under-reported disease due to failure to suspect, detect and diagnose the disease in many parts of the India. The Indian Council of Medical Research (ICMR) has been at the forefront of addressing public health challenges and in providing technical support to fill the aforementioned gaps in biomedical research. This Standard Operating Procedure (SOP) is a step towards such endeavor. It is a practical manual intended to be used as a reference guide by various laboratories in different parts of the country. The manual describes the safe and scientific methods of specimen collection, transport, culture, identification and molecular testing for diagnosing Anthrax.

I applaud the diligent efforts of the authors, editorial board members and the coordination team who have worked hard to bring this publication to its present form. I am hopeful that the users of this manual will be able to appreciate its value in their respective fields of work. The SOP will further evolve for clinical as well as research purposes through periodic revisions and updates.

I convey my best wishes to all.

Prof (Dr) Balram Bhargava

Secretary, Department of Health Research Director General, Indian Council of Medical Research (ICMR) Department of Health Research Ministry of Health and Family Welfare Government of India





भारतीय आयुर्विज्ञान अनुसंधान परिषद स्वाख्य अनुसंधान विमाग स्वाख्य एवं परिवार कल्याण मंत्रात्मय वी, वगदिगस्वामे भवन, अंसरी नगर नई दिल्ली - 110 029 (मारले) Indian Council of Medical Research Department of Health Research Ministry of Health & Family Welfare V. Ramalingaswami Bhawan, Ansari Nagar New Delhi - 110 029 (INDIA)

MESSAGE

The devastating effects of COVID-19 pandemic globally have taught usseveral lessonson public health preparedness to prevent, protect against and respond to health emergencies. Anthrax being a lethalbacterial disease with additional potential of Bacillus anthracis being used as a bioweapon, merits due attention from such perspective. It is also important to recognize that persistent hot spots of anthrax with repeated outbreaks have been reported in the pastfrom different parts of India. Awareness among the communities andtraining of healthcare workers and veterinary officersare therefore recommended as an essential component of preparedness for future emergencies.

Effective control of anthrax depends on timely diagnosis and reporting, which presently remains as major lacunae in several endemic regions of the country. This resource material has been developed as a reference manualto respond to such observations and to guide laboratories on standardized uniform protocolsfor diagnosis.

The team at the Indian Council of Medical Research (ICMR) has prepared this Standard Operating Procedure (SOP) in partnership with the relevant stakeholders. We believe that this SOP will be of tremendous use as a national guideline for handling suspected human, animal and environmental specimens towards diagnosing and managing anthrax. It will also add fillip to the existing biosafety measures in the country and facilitate responsible handling of agents of biothreat potential.

Dr Samiran Panda

Additional Director General and Head Division of Epidemiology and Communicable Diseases (ECD) Indian Council of Medical Research (ICMR) Department of Health Research (DHR) Ministry of Health and Family Welfare Government of India

Chief Patron

Dr Mansukh Mandaviya Union Minister of Health & Family Welfare

Patron

Dr Bharati Pravin Pawar Minister of State for Health & Family Welfare

National Advisors

- 1. Prof. Dr Balram Bhargava, Secretary, DHR and Director General, ICMR
- 2. Dr Samiran Panda, Additional DG and Head, ECD, ICMR

CONTRIBUTORS

Experts

- 1. Dr DCS Reddy, Ex Professor & HOD, PSM, IMS, BHU &Retd NPO, WHO-India
- 2. Dr Balaji Veeraraghavan, Professor, Department of Microbiology, CMC, Vellore
- 3. Dr John Antony Jude Prakash, Professor and Head, Department of Microbiology, CMC, Vellore
- 4. Dr MonilSinghai, Joint Director, Centre for Arboviral and Zoonotic Diseases (CAZD), National Center for Disease Control (NCDC), New Delhi
- 5. Dr Simmi Tiwari, Joint Director (PH), Program Officer, Division of Zoonotic Diseases programmes, National Center for Disease Control (NCDC), New Delhi
- 6. Dr Tushar Nale, Deputy Director, National Center for Disease Control (NCDC), New Delhi
- 7. Dr DivyaElangovan, Aravind Eye Hospital, Chennai
- 8. Dr SB Barbuddhe, Director, ICAR National Centre on Meat, Hyderabad, Telangana
- 9. Dr B.R. Shome, Director (Acting), ICAR-NIVEDI, Benguluru
- 10. Dr Satish B Shivachandra, Principal Scientist, ICAR-NIVEDI, Benguluru
- 11. Dr Mudassar Chanda, Principal Scientist, ICAR-NIVEDI, Bengaluru
- 12. Dr Deepak Rawool, Senior Scientist, ICAR National Research Centre on Meat, Hyderabad
- 13. Ms Saaya, Former Instructor, Department of Microbiology, CMC, Vellore

Editorial Committee

- 1. Dr Nabendu Chatterjee, Head, Division of BMS, ICMR
- 2. Dr Madhumathi J, Scientist C, Division of ECD, ICMR

ICMR Coordination Team

- 1. Dr Madhumathi J, Scientist C, Division of ECD, ICMR
- 2. Dr Lokesh Kumar Sharma, Scientist E, Division of ECD, ICMR

Administrative support

- 1. Ms Vandana Sundariyal, Technical Officer, Division of ECD, ICMR
- 2. Ms Tanvi Saxena, Data Entry Operator, Division of ECD, ICMR



CONTENTS

Minis

ry of H

alth & Family Welfare ment of India

			icma
Sl. No	Title	Page no.	INDIAN COUNCIL OF MEDICAL RESEARCH Serving the rulture sizes 121
I Uumon Anthro	N.V.		and week
1. IIuiiiaii Aiitiiia 1. Standard nr	LA acquitions	1	
1. Stanuaru pr 1 1	Safety measures	1	A A A A A A A A A A A A A A A A A A A
1.1	Storage instructions	1	
1.2	Disinfection	2	
1.4	Decontamination	$\frac{1}{2}$	
1.5	Spillage management	$\frac{1}{2}$	
1.6	Final disposal	3	
2. Laboratory	Diagnosis	4	
2.1	Specimen collection	4	
2.2	Specimen transport	5	
2.3	Microscopy	5	
2.4	Culture characteristics	7	
2.5	Biochemical Test	8	
3. Serology		8	
3.1	Direct immunofluorescence	8	
3.2	Detection of anti-PA antibodies	9	
4. Antimicrobi	al susceptibility testing	10	
4.1	Broth microdilution test	10	
4.2 5 Confirmatio	E-test method	12	
5. Commination	ol ol <i>Daculus aninracis</i>	14	
0. Anumerou 7 Molecular d	iagnosis	15	
7. Molecular u 7 1	Real time PCR for anthray	10	
7.1 7.2	Whole genome sequencing and analysis	10	
II Animal Anthra	v note genome sequeneing and analysis	17	
8 Introduction		20	
8 1	Specimen collection	20	
8.2	Specimen transport	20	
8.3	Specimen processing	$\frac{20}{20}$	
III Environments	al samnling		
9. Specimen co	llection	21	
9.1	Specimen processing	21	
9.2	Culture	21	
10. Appendix -	STAINS	22	
10.1	Gram stain	22	
10.2	Loeffler's Methylene blue	23	
10.3	Polychrome Methylene blue	23	
10.4	Malachite green staining	23	
11. Culture Me	edia	24	
11.1	Nutrient broth	24	
11.2	Nutrient agar	24	
11.3	Blood agar	24	
10 D' 11.4	Muller Hinton blood agar	24	
12. Biochemica	us Conhohydroto forme antation toot	25	
12.1	Carbonyurate termentation test	25	
12.2	Calatinase test	25	
12.5	Nitrate reduction test	20	
13. Reference		28	

13. Reference

ICMR Standard Operating Procedure for Anthrax

I. Human Anthrax

According to the Centers for Disease Control and Prevention (CDC) *Bacillus anthracis* Classified as a category Abioterrorism agent.

1. Important standard precautions

1.1. Safety measures :

All laboratory personnel should be well trained in biosafety and biosecurity before commencing work.

Laboratory personnel should ensure adherence to safety protocols

All laboratory procedures should be performed in the Biosafety Level 3 containment facility (BSL III) inside Biosafety cabinet Class III to minimize production of potentia aerosols.

General Good Laboratory Practices (GLP) must be followed as per the World Health Organization (WHO) guideline for Anthrax, 4th edition 2008.

Laboratory personnel should wear gown or laboratory coat with long sleeves and elastic cuffings and disposable gloves.

Appropriate disinfectant (usually 10% hypochlorite solution) should be prepared freshly on regular basis.

Screw capped non-breakable containers should be used for transport of specimens, cultures, etc.Autoclavablecarriers or secondary containers should be used for moving cultures within the laboratory.

1.2. Storage instructions :

Storage of samples and *B.anthracis* isolates in the laboratory are strictly prohibited

- Samples and isolates should be kept only till the results are reported.
- After reporting, the samples and isolates should be autoclaved and decontaminated as per protocol of the World Health Organization (WHO).
- Isolates may be sent to the reference laboratory for confirmations (National reference laboratory, Department of Clinical Microbiology, Christian Medical College, Vellore (CMC), Tamilnadu and Center forArboviral and Zoonotic Diseases (CAZD), National Center for Disease Control (NCDC), New Delhi).





1.3. Disinfection

- Laboratory accessories like pipettes, tips, loops, spreaders etc. should be autoclaved (e.g. in strong autoclavable bags) or fumigated or should be fully immersed in disinfectant (10% hypochlorite or formalin).
- Biosafety cabinet should be disinfected after use with 10% sodium hypochlorite solution.
- Contaminated itemsshould be kept in strong leak proof containers preferably within autoclavable bags.
- Hands should be thoroughly washed with soap and water and driedbefore leaving the facility.

1.4. Decontamination :

- All materials used including labware should be decontaminated by autoclaving at 121°C for 30 minutes to 1 hour preferably followed by incineration.
- Microscopic slides, cover slips and other sharp items should be placed in autoclavable sharp containers and autoclaved, preferably followed by incineration.
- Infectious disposable waste should be autoclaved followed by incineration.
- Non-autoclavable materials should be disinfected or fumigated.
- Laboratory clothing should be autoclaved before being sent to the laundry.

1.5. Spillage management :

- If accidental spillage or exposure to anthrax sporesoccurs in the laboratory (including spills of infectious substances like blood or body fluids), PPE should be worn before cleaning the spillage (Gloves, protecting clothing including face and eye shield).
- The spill should be covered with absorbent cloth or paper towels soaked with disinfectant (10% hypochlorite).
- The disinfectant should be applied concentrically beginning at the outer margin of the spill area and working towards the centre.
- The spill should be allowed to soak in the disinfectant for about one hour.
- The materials should be cleared using forceps. Broken glass or other sharps should be placed in autoclavable sharps container for disposal.
- After autoclaving, the material is incinerated.



1.6. Final disposal

- All cultures of *Bacillus anthracis*, specimens and disposable labware used for isolation, identification and performing molecular characterization should be decontaminated by placing in double autoclave bags for autoclaving followed by incineration.
- The disposable material should be autoclaved twice followed by incineration before final disposal.
- Fumigation should be carried out for non-disposable items.
- The decontaminated material should be sent for final disposal after securely bagging the contents kept in a covered trolley to the disposal area.

Instruction for disposal of Anthrax Carcasses (Refer Guidance Document Repository (GDR), Canada)

- AVOID USING LIME or other calcium products on carcasses or contaminated ground.
- Carcasses disposal should be done by incineration. However, deep burial is also an acceptable method.
- Performing an autopsy is prohibited when anthrax is being considered.
- Ensure all body openings (e.g. anus, mouth, nose etc.,) are plugged with absorbent material, non-perforated towel, cloths etc., to prevent leakage.
- Head should be covered with heavy duty plastic bag and tied to prevent scavenging and spreading of spores by insects, birds and animals.
- Finally, incineration should be done.
- Carcasses inadvertently open for post mortem examination or scavenging, it can be disposed by burning C Pyre or Pit/Burn Pits.
- Ensure that an adequate amount of fuel is added to completely reduce the carcass to ash.
- If incineration is not possible deep burial may be a viable option.

Incase of prolonged rains or logistical problem, carcasses disposal may be delayed. Under these circumstances, cover the carcass and surround the area with 10% formalin or hypochlorite solution. This should be protected from scavenging and then finally incinerated.



Note to the Clinicians :

Clinicians are encouraged to provide adequate volume of specimen and a second specimen obtained to improve diagnostic yield.

Swabs are generally discouraged. When adequate purulent material or vesicular fluid is available, as they carry very less material. Two swabs are needed, one for smear and the other for culture.

Clinicians are encouraged to contact the lab prior to sending suspected anthrax specimens, so that the lab persons take appropriate protective precautionary measures and also to be ready to process sample the moment it is received.

2. Lab Diagnosis

Procedure:

2.1. Specimen collection :

Specimen: Pus (exudate) aspirate, swab, fluid from bleb, ascitic fluid, Peritoneal fluid, CSF, blood, sputum, food(Intestinal anthrax) and stool.

Specimen	Container	Expected sample amount *	
Aspirated pus (exudate)	Universal disposable sterile container (wide mouthed, screw capped, leak proof)	<u>≥</u> 0.5 ml	
Pus swab# (when pus cannot be aspirated)	Sterile tube with two sterile swabs	Two swabs; Each should be insinuated beneath the scab and rotated 2-3 times each	
Fluids (bleb, ascitic, peritoneal etc.)	Universal disposable sterile container (wide mouthed, screw capped, leak proof)	≥0.5 ml	
CSF	Sterile Eppendorf (given in the lumbar puncture set)	≥0.5 ml	< T
Blood	BacT/Alert blood culture bottles (FAN)	Maximum:10 ml	
Sputum & Gastric aspirate	Universal disposable sterile container (wide mouthed, screw capped, leak proof)	≥0.5 ml	
Stool	Universal disposable sterile container (wide mouthed, screw capped, leak proof)	Maximum: 1/3 of container	



2.2. Specimen transport :

All specimens are collected in appropriate leak proof containers and secured, the sample container is put in a clear plastic cover with biohazard label, sealed and immediately transported to the laboratory (hand carried) and handed over to the lab staff.

Processing of sample

The sample is processed in the bio-safety cabinet BSLIII type 2B2. Smears are prepared and cultures inoculated here, including the culture follow-up.

Personal Protective Equipment (PPE) : Wear gloves, N95 mask, apron, visor and boot cover; The mandated Personal Protective Equipment (PPE) is to be worn during sample processing and follow-up.

2.3. Microscopy :

Direct microscopy :

Make 2 smears each on 3 separate glass slides, Perform Gram stain, polychrome methylene blue and methylene blue staining.

Gram stain :

Procedure :

Prepare the smear :

- 1. Use clean grease free and scratch free slides.
- 2. The slides are labeled with diamond marker or white pencil to ensure identification.
- 3. Use an inoculating loop or sterile pipette, or needle and syringe to place the liquid specimen (aspirate/fluid/pus/broth) at two point of a labeled slide.
- 4. Smear using swabs: Press the swab a few times onto the glass slide in a circular manner (impression smear on the surface of the slide).

Note: For specimens from suspected anthrax cases, two smears are prepared/clean glass slide.

- 5. To prepare smears from colonies growing on agar media by placing a small drop of saline onto the center of a labeled slide. Touch the center of a colony with a sterile smear loop and transfer a small amount of the bacterial colony to the drop of saline and mix the bacteria with the saline. Spread the mixture over an area approximately 2 sq.cm.
- 6. Allow smear to air dry.
- 7. Heat fix the smear. Mark the position of the smear as the reverse of the slide.
- 8. Note: If lot of purulent material is present, use methanol to fix the smear
- 9. Flood the fixed smear with the crystal violet allowing the stain to remain for one minute.
- 10. Pour off crystal violet and rinse slide in water using a wash bottle or gently flowing tap water until clear.
- 11. Flood the slide with iodine solution for one minute; rinse with gently flowing tap water or water from a wash bottle.
- 12. While holding the slide at a tilted angle, apply a few drop of decolorizer to the upper end of the slide, and allow the decolorizer solution to flow over the smear for 2 to 5 seconds.





1. Stop the decolorization after 2 - 5 seconds with a gentle flow of water. Do not apply decolorizer until the color stops running off because that will over decolorizer the bacteria.

- 2. Rinse with gently flowing tap water or water from the wash bottle.
- 3. Apply the counter stain safranine for 30 seconds and remove the excess stain by washing with gentle flowing tap water or water from the wash bottle water from a wash bottle.
- 4. Drain slide or air dry do not blot dry.
- 5. Examine the smear microscopically using the oil immersion objective (100x).

Gram stain report:

- Gram positive large thick rectangular bacilli
- Size 4x1µm
- Box car Shaped
- The ends of the bacilli are truncated giving a "Bamboo-stick" appearance
- Spores-Non bulging oval central or sub terminal spores
- In tissue single or short chains of bacilli without spores
- In culture long chains spores are seen

Methylene Blue Staining :

Purpose:

Methylene blue is a simple direct stain used to reveal the morphology of bacteria and to visualize the presence of capsule.

Expected Result:

Bacteria stain deep blue and leukocytes stain blue against a light blue staining background. Capsule of *B. anthracis* stains pink (this may not be seen often)

Polychrome Methylene blue (PCMB) staining protocol : Purpose :

Demonstration of capsulated *Bacillus anthracis* in clinical specimens from patients suspected to have anthrax or after capsule induction test.

Mc Fadyean's reaction :

In Polychrome methylene blue stains *B. anthracis* capsule appears as amorphous purple material surrounding blue bacilli known as (Mc Fadyean'sreaction).

Media Inoculation and Culture Characteristic : Media inoculation :

Nutrient agar (NA), Blood agar (BA), Nutrient broth (NB) (seal all the culture plates and tubes with adhesive tape).

Incubation : Blood agar (BA) 5-10% CO₂ incubator at 37°C Nutrient agar (NA), Nutrient broth (NB) aerobic incubation at 37°C

Duration of incubation :

Plates (BA & NA): 48 hours









Government of India

2.4. Culture characteristics : Nutrient agar

Colonies are 2-3mm in size irregular, raised dull grevish whitewith frosted glass (ground glass) appearance.

Medusa head appearance :

When colonies are viewed low power microscope the edge of the colony which is composed of long interlacing chains of bacilli appears as locks of matted hair, giving medusa head appearance, may have fringed edge or put out curled protrusions (tailing).

Colonies have tacky consistency.

Blood Agar:

Non-hemolytic dry white or greyish white wrinkled irregular, low convex colonies with tailing along the streak line.

Blood Agar with penicillin : Concentration of penicillin 0.05-0.50 units

Solid medium with Penicillin: Colonies have a string pearl appearance look due to the cells becoming larger and spherical because of the weaker cell walls under the action of penicillin, and cells tend to occurs in chain on surface of agar.

Nutrient broth : Floccular deposit with supernatant clear or slightly turbid

Do a hanging drop preparation from the overnight incubated nutrient broth (NB).

B anthracis is non-motile.

Note : *B* anthracis will be often observed in chains in NB, look for motility amongst bacilli which are singly placed or in pairs before reporting. If in doubt, re-inoculate into a fresh NB, incubate for 2 hours and then do a hanging drop for motility. Keep the slide in a petridish and autoclave (DO NOT discard the cover-slip in discard jar).

Capsule induction test :

The capsule induction test is performed for quick confirmation of B. anthracis from Blood agar (BA) & Nutrient agar (NA) plates.

When is it performed: When long, thick Gram Positive Bacilli with sharp ends, with spores, are seen in Gram stained culture smears.

Procedure :

Inoculate a single colony into 5 ml of sterile sheep blood, seal with adhesive tape and incubate in 5% CO₂ incubator for 12-24 hrs.

- 1. After 12-24 hrs incubation make smear of the blood using swab.
- 2. Discard the swab in 10% hypochlorite solution but NOT IN LYSOL.
- 3. Stain with polychrome methylene blue & methylene blue stain and look for the McFadyean's reaction.

Note: After the completion of the test autoclave the test tube with blood.

Susceptability to penicillin G:

- 1. Make a lawn culture of *Bacillus anthracis* on nutrient agar or blood agar.
- 2. Place 10 units of penicillin disk and incubate at 37°C for 16-18 hrs.
- 3. Zone of inhibition will be visible around the penicillin disk confirms susceptibility to penicillin.













Culture Smear :

Gram-staining : Reveals bamboo stick appearance, i.e. long chain of gram-positive bacilli with non-bulging spores (appears as empty space)

Spores:





They can be demonstrated using special stains, such as hot 5% malachite green(Ashby's method) or 0.25% sulphuric acid used in acid fast staining for spores. Spores appears as pink in colour bacilli appears as blue.

2.5. Biochemical Test :

Catalase Test :

- Catalase produce should be done inside the petri plate **Gelatin Liquefaction Test:**
- Slow liquefaction of gelatin Growth occurs as "Inverted fir tree appearance"(due to liquefaction of gelatin which occurs maximum at the surface, and then slows down towards the bottom)

Nitrate Test:

• Nitrate are reduced to nitrite.

Carbohydrate fermentation Test:

- Glucose, Maltose, Sucrose, Trehalose dextrin fermented without gas.
- Salicin not fermentated.

Safety Measures :

- Place petri dishes (or other culture containers) in a purpose-designed carrier or secondary container, such as a sandwich box, for movement around the laboratory. The carrier or container should be labelled with the agent, the operator's ID and date.
- Discard the plates/tubes into autoclave bags. Autoclave, preferably followed by incineration.
- Discard used slides in 10% hypochlorite and autoclaved
- Other sharp items into the sharps container which is autoclaved and then, preferably, incinerated also.
- Incinerate/autoclave other used disposable items of equipment.
- Double Autoclave the recyclable item.
- Fumigate or otherwise decontaminate non-disposable items of equipment which cannot be autoclaved.
- Decontaminate the safety cabinet after use with 10% sodium hypochlorite solution appropriate disinfectants.

3. Serology :

Antibodies appear in convalescent sera and can be detected by ELISA or immunodiffusion in gel method.

Direct Demonstration :

3.1. Direct Immunofluorescence test :

If detect capsular and cell wall polypeptide antigens by using fluorescent tagged monoclonal antibodies. It is used for confirmation of the diagnosis during bioterrorism outbreaks.



3.2. Detection of anti-PA antibodies

Purpose:

To detect anti-PA (protective antigen) antibodies, which are suggestive of exposure to *Bacillus anthracis*

Principle:





Indirect ELISA; IgG antibody to protective antigen (anti-PA) if present in the serum will bind to the antigen coated on the micro-titre plate. The amount of antibody present will be proportional to the intensity of color generated when enzyme linked conjugate (anti-species antibody) is added.

Sample: Canine serum

Procedure: As per the prescribed kit insert

Kit used: Anthrax Protective Antigen IgG ELISA (Cat No.SE120147); Sigma-Aldrich, St Louis, MO, USA)

Calculations :

- 1. Check Calibrator Factor (CF) value on the calibrator bottle. This value might vary from lot to lot. Make sure you check the value on every kit.
- 2. Calculate the cut-off value: Calibrator OD x Calibrator Factor (CF).
- 3. Calculate the Ab (Antibody) Index of each determination by dividing the O.D. value of each sample by cut-off value.

Validation criteria: All 3 criteria given below have to be satisfied for every ELISA run before declaring test results.

The three validation criteria to be met are:

- 1. The O.D. of the Calibrator should be >0.250.
- 2. The Ab index for Negative control should be <0.9.
- 3. The Ab Index for Positive control should be >1.2.

Interpretation:

- <0.9 No detectable antibody to PA IgG by ELISA
- 0.91.1 Borderline positive. Follow-up testing is recommended if clinically indicated.
- >1.1 Indicative of vaccination, current or previous Anthrax infection.

Quality assurance:

- Kit controls (PC & NC) should give expected results
- Results of the in-house QC and split sample testing should be within range.

PPE : As serum is not containing anthrax bacilli (only antibodies), standard PPE should be worn while testing

4. Antimicrobial Susceptibility Testing : (AST)

Antimicrobial Susceptibility Test is performed and interpreted as per the protocol described in Table 2K of CLSI M-100 S-20 (Vol30; No.1); January 2010.

The gold standard for determining antimicrobial susceptibility for Anthrax is the conventional broth microdilution (BMD) and E-test method which is based on the Clinical and Laboratory Standards Institute guidelines.

(Note : DD (disk diffusion) is not the recommended method for *B. anthracis*)

MIC is performed for Penicillin, Tetracycline, Doxycycline, Ciprofloxacin and Levofloxacin. by broth microdilution (BMD) and E-test method

4.1. Broth microdilution test :

Broth microdilution tests are done in microtiter plates with a final volume of 0.1 ml. Requirements:

- 1. U-bottom or V-bottom 96 well microtiter plates
- 2. Sterile graduated pipettes 10 ml
- 3. Sterile Pasteur pipettes
- 4. Overnight culture of test and control organisms
- 5. Required antibiotic in powder form (either from the manufacturer or standard laboratory accompanied by a statement of its activity in mg/unit or per ml. Clinical preparations should not be used for reference technique)
- 6. Required solvent for the antibiotic
- 7. Sterile distilled water 500 ml
- 8. Suitable nutrient broth medium (cation adjusted Muller Hinton Broth Medium).
- 9. Micropipettes $20 200 \,\mu$ l, $10 \,\mu$ l with sterile tips

Procedure:

- 1. Prepare antimicrobial stock and serial two-fold dilution of the antimicrobial range to be tested in broth as in macrobroth dilution.
- 2. Add 0.1 ml of the antibiotic into the wells
- 3. Include one sterility well (only broth, no antimicrobial/inoculum) and one growth control well (broth with inoculum, no antimicrobial) in each tray
- 4. Inoculum preparation:
 a. Prepare 0.5 McFarland matched suspension of the organism either by direct suspension or growth method
 b. Dilute 1:20 in water-tween 80 diluent (2ml in 40 ml diluent)
- 5. Add 0.01 ml of the standardized inoculum preparation in all the wells except the sterility well using prongs or with micropipettes.
- 6. Seal the plate with a plastic tape or in a plastic bag
- 7. Check purity of inoculum by subculturing onto non-selective agar plate
- 8. Inoculum verification is done as follows:
 - a. Take 0.01 ml of the growth control and dilute in 10 ml sterile saline
 - b. Plate 0.1ml onto non-selective plate
 - c. Count colonies after overnight incubation
 - d. Approximate count of 50 indicates an inoculum density of 5 x 105 CFU/ml
- 9. Incubate the trays stacked no higher than four plates at 35 ± 2 oC for 16 to 20 hours in ambient air





Reading of plates:

- 1. Check purity of inoculum
- 2. Check sterility control well. It should be clear/non-turbid
- 3. Check growth control for adequate growth of at least 2mm button formation
- 4. The MIC is the lowest concentration of antimicrobial agent that completely inhibits growth of the organism in the microdilution wells as detected by the unaided eye
- 5. When single skipped well is seen, read the highest MIC
- 6. Do not report results if more than one skipped well is present
- 7. For trimethoprime and sulfonamides, read the end point at the concentration in which there is = 80% reduction in growth as compared to the control

Recording template for microbroth dilution MIC:







4.2. E-test method :

The E test (AB Biodisk) which is a quantitative method for antimicrobial susceptibility testing applies both the dilution of antibiotic and diffusion of antibiotic into the medium. A predefined stable antimicrobial gradient is present on a thin inert carrier strip. E test can be used to determine Minimum Inhibitory Concentration.

When this E test strip is applied onto an inoculated agar plate, there is an immediate release of the drug. Following incubation, a symmetrical inhibition ellipse is produced. The intersection of the inhibitory zone edge and the calibrated carrier strip indicates the MIC value over a wide concentration range (>10 dilutions) with inherent precision and accuracy.

Requirements:

- 1. E-test strips
- 2. Media (Mueller Hinton Blood Agar)
- 3. 0.5 McFarland turbidity standard
- 4. Cotton swabs
- 5. Forceps/e-test applicator
- 6. Overnight culture of QC and test organisms

Procedure :

- 1. Remove e-test strips from freezer and allow to come to room temperature.
- 2. 0.5 McFarland standard matched suspension of the organism is prepared by growth method or direct colony suspension.
- 3. Streak a lawn culture of the standardized inoculum as in disk diffusion testing within 15 minutes of inoculum preparation
- 4. Allow excess moisture to be absorbed into the media
- 5. Apply strips onto agar surface with the forceps making sure the graduated surface is facing up.
- 6. Place only one strip for 90 mm petridish.
- 7. Do not change position of the strip once placed on the plate.
- 8. Incubate at $35\pm 2^{\circ}$ C for 16 to 20 hours in ambient atmosphere.



Fig Placement of multiple e-test strips on 150 mm and 100 mm petridishes



Reading:

- 1. Read MIC as the point where the inhibition ellipse intersects the scale.
- 2. Refer e-test reading guide from AB





E-test reading guide from AB Biodisk



Interpretative criteria:

Anti-microbial agent used	Susceptible MIC (µg/ml)	Non-susceptible MIC (µg/ml)	Comments
Penicillin	<u>≤</u> 0.12	<u>></u> 0.25	Can be extrapolated to amoxicillin
Tetracycline Doxycycline	≤1 ≤1		Non-susceptible strains are
Ciprofloxacin Levofloxacin	≤0.25 ≤0.25		very rare*

*If MIC value suggests isolate is not susceptible, the following has to be done

- 1. Confirm the identity of the isolate (as *B. anthracis*)
- 2. Confirm the AST results
- 3. Save the isolate
- 4. Submit the isolate to a reference laboratory for confirmation of ID & susceptibility

5. Confirmation of *B. anthracis*:

- 1. "Bamboo stick" appearance in Gram stain (GPB with non-bulging spores in culture isolates)
- 2. MacFadyean's reaction in polychrome methylene blue stain
- 3. Non-motile bacilli
- 4. Non-haemolytic colonies on blood agar after overnight incubation
- 5. Medusa head appearance on NA/BA
- 6. Capsule induction test positive
- 7. Penicillin susceptible (often)
- 8. PCR confirmation (detection of PX01 and PX02 plasmid)

6. Antimicrobial resistance in Anthrax

Rapid initiation of appropriate antibiotic therapy is crucial in the treatment of anthrax particularly in the systemic manifestations. Penicillin had been the drug of choice for all types of anthrax since 1940s, but naturally occurring strains are increasingly reported to be resistant. B. anthracis is sensitive to a broad range of antibiotics including tetracyclines, macrolides, aminoglycosides, fluoroquinolones, carbapenems, linezolid, clindamycin, rifampin, quinupristin-dalfopristin, daptomycin, and first-generation cephalosporins. However, second and third generation cephalosporin, cotrimaxozole to be avoided in the management of Anthrax, as they are relatively resistant to these drugs[1]

Most strains of naturally occurring *B. anthracis* have a chromosomally mediated, weak, inducible Beta-lactamase and cephalosporinase, and there have been rare reports of the development of resistance during therapy with penicillin, especially if subtherapeutic doses may have been administered. Being the most potent bioweapon the organism is expected to be ineffective to penicillin treatment due to the reason that the penicillin resistance can be induced in laboratory conditions. [2] The guidelines for bioterrorism-associated anthrax recommend use of fluoroquinolones, carbapenems, and doxycycline until resistance testing is available, as â-lactam resistance in such strains is presumed to be likely to be ruled out.

The gold standard for determining antimicrobial susceptibility for Anthrax is the conventional broth microdilution (BMD) method which is based on the Clinical and Laboratory Standards Institute guidelines. [3] This method is growth- dependent and requires an incubation period of 16 to 20 h for B. anthracis. Other methods commonly used for AST are agar dilution and diffusion-based assays such as the disc diffusion and the Etest. These alternative methods, although easier to handle, require similar incubation times since visible growth is required for the interpretation of the results. A comparison of standard broth micro-dilution and E-test agar diffusion on *B. anthracis* isolates, conducted by the CDC [4,11], found that there was no statistically significant difference between both methods for any of the antimicrobial agents tested, except for penicillin in which the E-test method was two-fold dilutions lower than the standard broth microdilution method

Reports of antimicrobial susceptibility profile of *B. anthracis* are scarce. This is probably due to the less number of human anthrax cases reported. Long-term antibiotic therapy, as would be administered for post exposure prophylaxis in anthrax, might induce antimicrobial resistance by the selection of resistant mutants. The possibility of an inducible beta-lactamase activity should be taken into consideration under clinical conditions where high numbers of organisms are to be expected. Though fluoroquinolones were highly active, strains that are resistant to fluoroquinolone has been isolated by Price et al invitro studies due to development of mutations in gyrA, parC and gyrB [5].Beta-Lactam-resistant strains have been attributed to the derepression of cephalosporinase. Doxycycline resistance was conferred on B. anthracis by transfection with a pBC16 plasmid carrying a tetracycline resistant gene, tet [6]

Lightfoot et al[7] evaluated nine antimicrobial agents with 33 epidemiologically distinct isolates by the agar dilution method. 90% of the strains were determined to be sensitive to penicillin, amoxicillin, gentamicin, streptomycin, erythromycin, tetracycline, and ciprofloxacin. Doganay et al [8] also used the agar dilution method for testing thirty antibiotics with 22 B. anthracis strains. They identified additional beta-lactams, aminoglycosides, clindamycin, vancomycin, and ofloxacin as antibiotics having activity.



Table 1: Antibiotic susceptibility profile of various Human Anthrax outbreaks reported globally

		N 6	1	6 1 1 1			
Human	Number	Year of	ar of Susceptibility Profile			Ref	
Outbreak Cases	of isolates	Isolation	Pencillin	Fluoroquinolones	Erythromycin and	Tetracycline	
					Clindamycin		
Turkey	22	1981- 1988	100%	100%	-	-	[8]
Ukraine	Case report	1992	0%	100%	100%	100%	[12]
France	96	1994- 2000	85%	100%	95-100%	100%	[9]
India	Case report	1997	0%	100%	100%	100%	[13]
Romania	21	2000- 2004	90%	100%	100%	100%	[10]
Florida	11	2001	95%	100%	98%	100%	[3]
Atlanta, Georgia	65	2004	96%	100%	100%	100%	[11]
Zambia	5	2011	100%	100%	100%	-	[14]
Kyrgyzsta n, Central Asia	138	2005- 2015	78%	>95%		>95%	[15]

7. Molecular Diagnosis:

PCR with specific primers can be used for further confirmation. Molecular Typing: It is used for epidemiological studies to trace the source of infection.

Various methods available are as follows :

- MLVA (Multiple locus variable number of tandem repeat analysis)
- AFLP (Amplified fragment length polymorphisms)

7.1. Real-time PCR for anthrax

Purpose : To detect DNA of *B anthracis* directly form samples and from cultures

Principle :

A segment of DNA unique to the pXO1 and pXO2 gene of *B anthracis* is detected by realtime PCR using a commercial TaqMan assay kit from Thermo Scientific

DNA extraction :

Samples : Pus, purulent tissue, aspirates, fluids and stool

Amount of sample :

About 25 mg of purulent tissue (which is hand-ground in a glass mortar and pestle) or 200 μ l of aspirate or fluid specimen or 200mg of stool.

Pre-treatment (to improve DNA yield): Add 100 μ l of Lysozyme (Sigma Aldrich); Pulse vortex for 15 seconds; Incubate at 37° C for one hour.

DNA from tissues and aspirates/fluids or cultures is extracted using the DNeasy Blood & Tissue Kit or Blood Mini kit (Qiagen) as per manufacturer's instructions after the pre-treatment procedure

The QIAamp DNA Stool (Qiagen) protocol is used to extract DNA from lysozyme treated stool.

The spin column protocol is used for extracting DNA from all samples. DNA purity is measured using 260/280 nm in a Thermo Scientific[™] NanoDrop 2000 and the ratio is expected to be above 1.5.

The eluted DNA is stored in 3 or 4 aliquots (20-25 μ l/ aliquot) at -70° C pending PCR amplification.

Real time PCR :

Kit based test used : TaqMan[™] Bacillus anthracis Detection Kit (Applied Biosystems; Part No. 4382486)

Targets amplified :

The targets amplified are pXO1 and pXO2 gene respectively and an internal positive control (IPC)

Type of assay : TaqMan

Components of the kit are as follows : Kit to be stored 20°C

- 1. Bacillus anthracis Target Assay Mix, pXO1: Purple cap; 300 µL (1 vial)
- 2. Bacillus anthracis Target Assay Mix (10X), pXO2: Green cap; 300 µL (1 vial)
- 3. Negative control: White cap; $1000 \ \mu L (1 \ vial)$
- 4. Fast PCR Master Mix (2X): Red cap; 750 μL (2 vials)

Note : pXO1 & pXO2 probes are labelled with FAM while the IPC is labelled with VIC





Real-time PCR equipment used :

- Real time PCR Machines
- All samples and controls are run in duplicate Workflow as described in the kit insert
- Prepare PCR master mix

Component	Contains	pXO1 master mix/reaction	pXO2 master mix/reaction
Fast PCR Master Mix (2X)	Polymerase & PCR buffer with MgCl ₂	15 µl	15 µl
Target Assay Mix (10X)	Forward & reverse primer plus TaqMan probe	3 µl	3 µl
Total volume		18 µl	18 µl

Note : When preparing the master-mix, always prepare for n+1 reactions. This will ensure adequate master mix is dispensed for all reactions.

1. Example if, 10 samples including controls have to be run including controls, master mix is prepared for 22 reactions (as all samples and controls are run in duplicates).

- Distribute 18 µl into each 0.1 ml tube or strip
- Set up amplification (create and store run file before-hand) PCR amplification parameters

	Stages	Temperature	Time	Cycles
Stage 1	Enzyme activation	95°C	20 seconds	1 cycle
Stage 2 (PCR)	Denaturation	95°C	3 seconds	45 gyales
	Annealing & extension	60° C	30 seconds	45 Cycles
Stage 3	Hold	4°C	×	

• Amplification performed on "Fast" run mode

• After the run is completed (sample temperature has reached 4° C) the amplification plots for all samples are viewed

2. Set the baseline: Baseline is set by default above 0.1

3. Check each sample for a FAM[™] dye (target-specific) signal and a VIC[®] dye (IPC) and note the Ct value

The threshold Ct value for a positive result is <35 for both duplicates. A negative result is given when the Ct value is >40 for both duplicates

If value is between 35 and 40 for any sample, assay is repeated, and if Ct<35 it is reported as positive. If Ct is still between 35 & 40 it is reported as borderline signal.





Interpretation

FAM dye signal (target)	VIC® dye signal (IPC)	Result	
+	+	Positive	
+	-	Positive	
-	+	Negative	
-	-	Troubleshoot*	

Note : * Please refer kit insert or manual

7.2. Whole genome sequencing and analysis

After DNA extraction as described in SOP 3, whole genome analysis will be performed on the (MinIon Nanopore Technologies) to obtain long sequence, fairly accurate reads. Sequencing will also be performed month (HiSeq Illumina) to obtain short but highly accurate reads. Assembly and analysis will be performed using the standard tools and pipelines used in the department for this purpose.

II. Animal Anthrax

8. Introduction

Zoonotic disease :

- Herbivorous animals such as cattle, sheep and Elephant less often horses and pigs
- Infection occurs in susceptible animals by ingestion of spores present in the soil
- Animal develops Fatal septicemia, localized cutaneous lesions and discharge large number of bacilli from mouth, nose and rectum
- These bacilli sporulate in soil and remains as source of infection for man

8.1. Specimen collection and Transport :

Legislation in most countries forbids postmortem examination of animals that have died of anthrax. Animals that have died suddenly and unexpectedly should not be necropsied unless anthrax has been ruled out as the cause of death

- Blood stained fluid may exude from one or more body orifices.
- Petechiae and ecchymoses are often present in unpigmented or hairless areas of the skin.
- Blood from veins from ear pinnae
- Swabs from the lesions

8.2. Specimen for transport ("double-bagging")

- The specimens should be collected into sterile containers The containers should be wiped down with hypochlorite (10 000 ppm) and, with outer gloves changed first, put into an outer, secondary container (double-bagged). If the secondary container is a plastic bag, then this should be of good quality. It should, in turn, be sealed and, for transport, be put into a good-quality cool box or a strong plastic or metal container with a lid that can be made secure.
- The secondary and outer containers should bear the relevant hazard labels.

Generally, specimens should be stored at 28 °C. Preferably they should be transported in cool boxes, especially in hot weather and when the time interval between collection and delivery to the laboratory is likely to be more than 12 hours.

8.3. Specimen Processing

The sample is processed in the bio-safety cabinet type 2B2 in the Mycobacteria Lab (BSL III). Smears are prepared and cultures inoculated here, including the culture follow-up.

Personal Protective Equipment (PPE): Wear gloves, N95 mask, apron, visor and boot cover; The mandated Personal Protective Equipment is to be worn during sample processing and follow-up

For microscopy culture and confirmation of anthrax refer Human Anthrax SOP (refer pg.)







III. Environmental sampling:

9. Sample collection :

9.1. Specimen Processing :

Sample collected : Soil sampling

Sample collection areas : The soil samples in anthrax endemic areas especially in pastures where cattle graze (in endemic/non-endemic/suspected anthrax is reported) will be sampled for anthrax using the more efficient ground anthrax bacillus refined isolation (GABRI) procedure, described by Fasanella et al.

When a suspected outbreak is reported, soil samples corresponding to the head and tail area of the carcass should be collected. These are spots most likely to have been contaminated by spore contaminated blood (Braun et al 2022). Using appropriate PPE, samples should be scooped from the surface (50-70g; ~ half full 50 ml conical tube/flask).

Data collected : Data on soil type, pH, precipitation and temperature should be noted inaddition to the GIS co-ordinates.

9.2. Culture :

Soil processing for culture :

Ground anthrax bacillus refined isolation (GABRI) procedure (as described by Fasanella et al, 2013):

- 1. A 7.5g aliquot of soil will be added to 22.5 ml of washing buffer consisting of deionized water containing 0.5% Tween 20.
- 2. After 30 minutes of washing by vortexing, the suspension will be centrifuged at 2000 rpm for 5 min to eliminate gross debris.
- 3. The harvested supernatant is incubated, aerobically, at 64°C for 20 min to eliminate vegetative forms of B. anthracis.
- 4. After incubation, 5 ml of supernatant is added to 5 ml of Tryptose Phosphate Broth containing 125 ig/ml of Fosfomycin.
- 5. Then, from each sample, 10 plates of TMSP are seeded with 1 ml/plate of the mix and incubated, aerobically, at 37°C.

6. **Medium Use :**

TSMP medium (Trimethoprim Sulfamethoxazole Methanol (5 ml/lt) PolymyxinMedium): The TSMP is nothing but Columbia blood agar with trimethoprim (16 mg/lt), sulfamethoxazole (80 mg/lt), methanol (5 ml/lt) and polymyxin (300,000 units/lt).

PLET medium It consists of polymyxin, Lysozyme, Ethylene diamine tetra acetic acid (EDTA) and Thallons acctate added in heart infusion agar. It has been devised to isolate *B.anthracis* from mixture of other spore forming bacilli.

- 7. After 24 and 48 hours of incubation, each plate is examined and the colonies of *B.anthracis* appears 2-3 mm rough, circular creamy white with ground glass appearance.
- 8. Colony morphology of *B. anthracis* should be follow as per the SOP in human anthrax





10. Appendix

10.1. Gram Stain :

Purpose:

Gram stain is most commonly used differential stain which divides bacteria into two major groups as Gram positive and Gram negative.

Ingredients and preparation :

Crystal Violet Stain:

Crystal violet	1.0 gm
Sodium bicarbonate 5%	1.0 ml
Distilled water	99.0 m

- Add 1gm of crystal violet into a mortar.
- Using a pestle grind it well.
- Then add sodium carbonate little by little to get a smooth paste.
- Finally add water and mix well.
- Filter through a filter paper into a stoppered bottle.

Gram's Iodine :

Iodine crystal	2.0 gm
1N Sodium hydroxide	10.0 ml
Distilled water	90.0 ml

- Add NaOH to the iodine crystals, kept in a mortar.
- Grind the paste to get a smooth paste.
- Add distilled water and mix well
- Filter through a filter paper into a stoppered brown bottle.
- Keep away from sunlight.

Acetone 100%

Safranine Stain :

Safranine	0.34 gm
Absolute alcohol/rectified spirit	10.0 ml
Distilled water	90.0 ml

- Add absolute alcohol to the powder kept in a mortar.
- Grind the paste into a smooth paste.
- Add distilled water and mix well
- Filter through filter paper into a stoppered bottle.

Use:

Widely used in diagnostic bacteriology to differentiate Gram positive and Gram negative organisms.

Quality Control :

Positive control : a heat fixed smear of *S. aureus*, Gram positive violet/purple cocci in clusters.

Negative control : a heat fixed smear of *E. coli*, Gram negative pink/red bacilli. If the control slide is improperly stained, repeat with a new control slide and a slide prepared



from clinical material. Because most problems are due to personal techniques, adjust staining time and technique. However, if proper results are not obtained, after those adjustments, consider a possible problem with the stains.



Notify supervisor if unable to obtain correct results.

10.2. Loeffler's Methylene Blue :

Ingredients :

Methylene blue	0.2 gm
Absolute alcohol/rectified spirit	10.0 ml
Distilled water	90.0 ml

Preparation :

- Weigh methylene blue and put it in the mortar
- Add absolute alcohol little by little and grind the powder into a smooth paste with the mortar.
- Add distilled water and mix it thoroughly.
- Filter through a filter paper into a stoppered bottle.

Use:

This simple stain is used to make out clearly the morphology of organisms e.g. *H.influenzae* in CSF, *N.gonorrhoeae* in urethral pus.

Quality Control :

Usea heat fixed smear of Staphylococcus aureus and H.influenzae

10.3. Polychrome Methylene Blue :

Preparation :

This is made by allowing Loeffler's methylene to? ripen' slowly. The stain is kept in bottles, which are half filled and shaken at intervals to aerate the contents. The slow oxidation of the methylene blue forms a violet compound that gives the stain its polychrome properties.

NOTE : The ripening takes 12 months, or more to complete or it may be ripened quickly by the addition of 1% potassium carbonate to the stain.

Expected Results :

The capsule of *B. anthracis* is seen clearly as pink amorphous material surrounding the blue - black bacilli (M'Fadyean reaction)

10.4. Malachite Green staining 5% :

- Make a smear and fixation can be done by heat or alcohol fix
- Place the slide over a beaker of boiling water, resting it on the rim with the bacterial smear uppermost;
- Cover with 5% aqueous solution of malachite green
- Stain for 5 minutes , adding more stain solution if the stain covering the smear starts to dry

Or

- Place the slide in a moist chamber (a petri dish with moistened filter paper will do)
- Cover with 5% aqueous solution of malachite green
- Leave to act for 60 minutes Then, following either procedure
- Wash of stain with water using wash bottle (into hypochlorite solution)
- Counterstain with 0.5% safranin or 0.05% carbolfuchsin for 30 seconds
- Wash again (into hypochlorite solution) and allow to dry
- Expected results : Spore appear green and the vegetative bacilli red

11. Culture Media :

11.1. Nutrient Broth (NB)

Peptone Beef extract (Lab Lemco) Sodium chloride Distilled water 1.0 gm 0.4 gm 0.5 gm 100.0 ml

Weigh out all the ingredients as above, peptone should be taken last, because it sticks to the paper on exposure. Mix the ingredients and dissolve them by heating. When cool, adjust the pH to 7.4 7.6. Distribute in tubes, bottles or flasks and sterilize by autoclaving.

Use:

This is a basal medium and is also used to grow non festidious organisms for various purpose.

11.2. Nutrient agar (NA)

1.5 to 1.8 gm
0.05 gm
0.05 gm
100.0 ml

Mix the agar in nutrient broth and heat to dissolve. When cool adjust the pH to $7.5 \ 0 \ 7.6$. Sterilize by autoclaving. Pour as plates or slopes. To make deeps, reduce agar concentration to 0.5%.

Use:

This is used as a base for many media. Only nonfastidious organisms will grow on this.

11.3. Blood agar (BA)

Sterile defibrinated sheep blood	10 ml
Nutrient agar (melted)	100.0ml

Pour about 10 ml of melted Nutrient agar, as a base, into sterile petri dishes and allow setting. This forms a thin base for pouring blood agar. Add steriledefibrinated sheep blood (5 10%) to nutrient agar, the latter should be cooled to about 45 50°C before blood is added. Mix well and pour about 15 ml of blood agar over the base in each petri dish. Human blood is not recommended for the preparation of blood agaras it may contain certain antibacterial substances and hence unsuitable for use in media preparation. Alternately blood agar may be made with no agar base.

Use:

It serves as an enriched medium and a differential medium for haemolytic organisms. Most common pathogens grow on it.

Quality control:

Viridans Streptococci- good growth with á-haemolytic colonies Group A *Beta-haemolytic Streptococci-* good growth with Beta-haemolytic colonies

11.4. Mueller Hinton Blood agar (MHA with 10% sheep blood)

MHA is made from commercially available dehydrated medium Difco TM MHA. Catalogue No: 225250, (BD, Sparks, MD, USA).

38 gms of the powder provided is dissolved in 1 liter of dist. water mixed thoroughly boiled for 1 minute to completely dissolve. Adjust pH to 7.3+0.2. Autoclave at 121° C at 15 mts.Do not over heat.



Approximate formula per liter is as given below.

2.0 gm
17.5 gm
1.5 gm
17.0 gm
1000 ml

Use:

Standard medium for antimicrobial susceptibility testing - fastidious organism.

12. Biochemicals

12.1. Carbohydrate Fermentation Test :

Purpose:

To determine the ability of an organism to ferment a specific carbohydrate incorporated in the basal medium, producing acid or acid and gas.

Test medium and materials :

Prepare sugar solutions as described below for different groups of organisms and dispense in 3 4 ml quantities into test tubes (12 x 100 mm). Introduce Durham's tubes into glucose broth for the detection of gas production. Autoclave at115°C for 10 minutes. Disaccharides like lactose and sucrose are better filtered and added to sterile basal medium. (Glucose, sucrose, maltose, trehalose, dextrin, and salicin)

, sucrose, manose, irenatose, uez	and same
Sugar	0.5 gm
Nutrient broth base	100.0 ml
Bromthymol blue indicator	
(0.2% alcoholic)	1.2ml

Procedure : Inoculate the test organism to the carbohydrate medium incubate at 37° C aerobically for 24-28 hours

Carbohydrate sugars	Positive control	Negative control
Glucose	ATCC 25922 E. coli	Acinetobacter baumannii complex
Sucrose	ATCC 25922 E. coli	Acinetobacter baumannii complex
Maltose	ATCC 25922 E. coli	Acinetobacter baumannii complex
Trehalose	ATCC 25922 E. coli	Acinetobacter baumannii complex
Dextrin	Cornybacterium diptheriae	Diphtheroids
Salicin	ATCC 700603 Kleb. pneumoniae	Acinetobacter baumannii complex

Result : Positive-acid production with yellow colour **Negative-** no change in colour

12.2. Catalase Test :

Purpose:

To detect the ability of an organism to produce the enzyme catalase which breaks down Hydrogen peroxide to water and nascent oxygen.

To differentiate between *Staphylococci* and *Streptococci*.



Test material :

- Growth on Nutrient agar.
- 3% hydrogen peroxide
- A young culture of *S. aureus* on NA agar.

Test procedure :

- Emulsify part of a colony in sterile saline on a clean glass slide.
- Add a drop of 3% H2O2.
- Inoculate a known positive control always

Reading and interpretation :

Appearance of gas bubbles immediately after the addition of H2O2 is to indicates the presence of catalase enzyme.

Eg: *S.aureus* positive *Streptococci* negative

Note:

- Donot use Nichrome wire as this may give false positive reaction
- Do not use growth on BA to test catalase as RBCs might give false positive reaction.

12.3. Gelatinase Test :

Purpose:

To test the ability of an organism to produce a proteolytic enzyme gelatinase which liquefies gelatin.

Test medium :

Gelatin medium :	
Nutrient broth	100.0 ml
Gelatin	6.0 gms
Adjust pH 7.2 7.4	-

- Heat ingredients other than gelatin, to dissolve in a water bath.
- Then add gelatin little by little to get a uniform solution.
- Dispense 5.0 ml per tube.

Autoclave at 121°C for 15 mts.

Test procedure :

- Inoculate the test organism as a stab into nutrient gelatin.
- Incubate at 37°C for 48hrs to a few days.

Reading and interpretation :

Keep the gelatin cultures in a beaker of ice or in the refrigerator before taking reading. Liquefication of gelatin i.e. absence of setting of gelatin at low temperature indicates gelatinase activity.

E.g., *S. aureus*, *Pseudomonas spp*. Positive *S. epidermidis*, *E.coli* Negative

Note:

Gelatin normally melts at about 24 - 27°C and sets below 20°C.



12.4. Nitrate Reduction Test :

Purpose:

To test the ability of an organism to reduce nitrate to nitrite and occasionally to gas, molecular nitrogen.

This is used for the identification of Enterobacteriaceae.

Test medium and reagents : Potassium Nitrate Broth :

0.2 gm
5.0 gm
1000 ml

Transfer into tubes in 5 ml amounts and autoclave.

Test reagent :

Solution A:	
Sulphanilic acid	8.0 gm
5N acetic acid	1000 ml
Dissolve sulphanilic acid in acetic acid.	

Solution B :

Beta-naphthylamine	5.0 gm
5N acetic acid	1000 ml
Dissolve Alpha-naphthylamine in acetic acid.	

Test procedure :

- Inoculate the organism into potassium nitrate medium.
- Incubate at 37°C for 24 96 hrs.
- Immediately before use, make the test reagent by mixing equal volumes of solutions A and B.
- Add 0.1 ml of the test reagent to the test culture.

Reading and interpretation :

A red color developing within a few minutes indicates the presence of nitrite and hence the ability of the organism to reduce nitrate to nitrite.

E.g. Enterobacteriaceae

If red color does not develop, a pinch of zinc powder is added and shaken. If a red color develops on addition of zinc, it confirms that the test was negative.

E.g. Erwinia spp

If no color develops following addition of Zinc, it confirms that nitrite was further reduced to nitrogen.

E.g. Ps. aeruginosa.



13. Reference :

Dr. Naresh Jindal .Principal Scientist, Veterinary Public Health & Epidemiology, LalaLajpatRai University of Veterinary & Animal Sciences (LUVAS). Hisar 125004, Haryana. <u>nareshjindal1@gmail.com</u>

Dr. S ByregowdaDirector, Institute of Animal Health and Veterinary BiologicalsHebbal, Bangaluru 560024. <u>smbyregowda.iahvb@gmail.com</u>

Dr. Raveendra Hedge. Associate Professor, IAH & VB, Hebbal, Bangaluru: ravihegde63@gmail.com

Dr. ShivarajMura. Assistant Professor, Department of Diagnostic Bacteriology and Mycology, SRDDL, IAHVB, Bangaluru. <u>shivaraj.murag@gmail.com</u>

Dr. NiranjanSahoo. Professor/Head, Epidemiology & Preventive Medicine cum Coordinator, Centre for Wildlife Health, Odhisa University of Agriculture & Technology (OUAT), Bhubaneswar, Odisha- 751003. <u>niranjansahoo@hotmail.com</u>

Dr. D. Rani Prameela. Professor & Head, Department of Veterinary Microbiology, State Level Diagnostic Laboratory, Sri Venkateswara Veterinary University, Tirupati 517502, Chittoor District, Andhra Pradesh. <u>raniprameela.dr@gmail.com</u>

Dr. M. VivekanandanPillai.Professor/Head, Department of Emergency medicine and Trauma, Jawaharlal Institute of Post graduate Medical Education & Research, Pondicherry. <u>drvivek@gmail.com</u>

Dr. SanghamitraPati. Scientist G and Director, Regional Medical Research Centre, Chandrasekharpur, Bhubaneswar, Odhisa, India. <u>sanghamitra.pati@icmr.gov.in</u>

Braun P, Beyer W, Hanczaruk M, Riehm JM, Antwerpen M, Otterbein C, Oesterheld J, Grass G. Reoccurring Bovine Anthrax in Germany on the Same Pasture after 12 Years. J ClinMicrobiol. 2022 Mar 16;60(3):e0229121. doi: 10.1128/jcm.02291-21. Epub 2022 Mar 16. PMID: 35195442; PMCID: PMC8925895.

Fasanella A, Di Taranto P, Garofolo G, Colao V, Marino L, Buonavoglia D, et al. Ground Anthrax Bacillus Refined Isolation (GABRI) method for analyzing environmental samples with low levels of Bacillus anthracis contamination. BMC Microbiol. 2013 Jul 18;13:167

Rotem S, Steinberger-Levy I, Israeli O, Zahavy E, Aloni-Grinstein R. Beating the Bio-Terror Threat with Rapid Antimicrobial Susceptibility Testing. Microorganisms. 2021 Jul;9(7):1535.

Dutta TK, Sujatha S, Sahoo RK. Anthrax-update on diagnosis and management. J Assoc Physicians India. 2011 Sep 1;59(1):573-8.

CLSI. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria, 3rd ed.;

CDC. Investigation of bioterrorism-related anthrax and interim guidelines for exposure management and antimicrobial therapy. MMWR 2001;50:909-919





Price, L. B., Volger, A., Pearson, T. et al. (2003). In vitro selection and characterization of Bacillus anthracis mutants with highlevel resistance to ciproflfloxacin. Antimicrobial Agents and Chemo therapy 47, 23625.

Chen, Y., Succi, J., Tenover, F. C. et al. (2003). b-Lactamase genes of the penicillinsusceptible Bacillus anthracis Sterne strain. Journal of Bacteriology 185, 82330

Lightfoot NF, Scott RJ, Turnbull PC. Antimicrobial susceptibility of Bacillus anthracis. Salisbury Med Bull Suppl 1990;68:95-98.

Doganay M, Aydin N. Antimicrobial susceptibility of Bacillus anthracis. Scandinavian journal of infectious diseases. 1991 Jan 1;23(3):333-5.

Cavallo JD, Ramisse F, Girardet M, Vaissaire J, Mock M, Hernandez E. Antibiotic susceptibilities of 96 isolates of Bacillus anthracis isolated in France between 1994 and 2000. Antimicrobial agents and chemotherapy. 2002 Jul;46(7):2307-9.

Caplan DM, Ivana S, Caplan ME. Susceptibility to antibiotics of Bacillus anthracis strains isolated in Romania. Roumanian archives of microbiology and immunology. 2009 Apr 1;68(2):106-10.

Mohammed MJ, Marston CK, Popovic T, Weyant RS, Tenover FC. Antimicrobial susceptibility testing of Bacillus anthracis: comparison of results obtained by using the National Committee for Clinical Laboratory Standards broth microdilution reference and Etest agar gradient diffusion methods. Journal of Clinical Microbiology. 2002 Jun;40(6):1902-7.

Bradaric N, Punda Polic V. Cutaneous anthrax due to penicillin resistant Bacillus anthracis transmitted by an insect bite. The Lancet 1992;340:306-7.

Lalitha MK, Kumar A. Anthrax- A continuing problem in South India. Indian Journal of Medical Microbiology 1996;4:63-72

Hang'ombe MB, Mwansa JC, Muwowo S, Mulenga P, Kapina M, Musenga E, Squarre
D, Mataa L, Thomas SY, Ogawa H, Sawa H. Human animal anthrax outbreak in the Luangwa valley of Zambia in 2011. Tropical doctor. 2012 Jul;42(3):136-9.

Kutmanova A, Doganay M, Zholdoshev S. Human anthrax in Kyrgyz Republic: Epidemiology and clinical features. Journal of Infection and Public Health. 2020 Aug 1;13(8):1161-5.